



TECHNICAL DATA SHEET

chemical indicator strips for steam sterilization Class 4



Steam indicator strips
Model No.: YSIS-15100
Size: 15mm × 100mm

Features:

- Be consisted of special card, standard color patch, and indicator strip
- To display indicator color patch (blue) by printing indicator ink onto a strip of special paper (standard color patch).
- Lead free & Latex free & non-toxic

Standards:

YIPAK chemical indicator strips is classified as Class 4 according to ISO 11140-1 standard. CE marking is not applicable for indicator products, because this product is not considered a medical device under the clause of the EC-directive 93/42/EEC.

Intended Use:

YIPAK Indicators are allow the end user to identify incorrect packaging, loading problems and sterilization failures

Operation Conditions:

121°C 15-20mins, 134°C 3.5-5mins, the initial color blue changed to black or lighter after a complete steam sterilization process

Storage:

Keep stored in the dark at room temperature 15°C~30°C and 50% relative humidity, avoid to contact with corrosive gases.

Shelf time:

Since the date of production 3 years

Packing:

YIPAK indicator strips are packed in recyclable, brown hard carton cases, consisting smaller units of 200 strips in each box.

Each case bears a label with the necessary information/instructions of the content in accordance with ISO 11140-1:2005.